



HEALTHCARE

In Healthcare, temperature management is crucial. Not only is food safety of the utmost importance, but so is preserving the efficacy and integrity of your vaccines, medications, blood, and test samples. This is quite the challenge, especially with shrinking budgets and constricting government regulations. SMART Temps® can help.

SMART Guard®

A wireless thermometer for automatic temperature monitoring of cold and hot equipment.



- Monitors cold/hot storage equipment 24/7/365.
- Identifies potential equipment maintenance issues.
- Sends alerts via email, text, or phone if temperatures are not in the proper range.
- Onboard memory eliminates potential lost or missed temperatures.

| Site | Step | Headend | Item | Leftover | Temp | Time | LBI | Result | Action | User |
|------|------|---------|------|----------|------|---------------------|-----|--------|--------------|------|
| ... | ... | ... | ... | ... | ... | 06/19/2012 11:47 am | ... | Pass | Keep Cooling | User |
| ... | ... | ... | ... | ... | ... | 06/19/2012 11:42 am | ... | 1 | Pass | User |
| ... | ... | ... | ... | ... | ... | 06/19/2012 11:41 am | ... | 1 | Pass | User |
| ... | ... | ... | ... | ... | ... | 06/19/2012 11:41 am | ... | 1 | Pass | User |
| ... | ... | ... | ... | ... | ... | 06/19/2012 11:36 am | ... | 1 | Pass | User |
| ... | ... | ... | ... | ... | ... | 06/19/2012 11:32 am | ... | High | Refrigerate | User |
| ... | ... | ... | ... | ... | ... | 06/19/2012 11:32 am | ... | 1 | Pass | User |
| ... | ... | ... | ... | ... | ... | 06/19/2012 11:30 am | ... | 1 | Pass | User |
| ... | ... | ... | ... | ... | ... | 06/19/2012 11:27 am | ... | 1 | Pass | User |
| ... | ... | ... | ... | ... | ... | 06/19/2012 11:27 am | ... | 1 | Pass | User |
| ... | ... | ... | ... | ... | ... | 06/19/2012 11:05 am | ... | 1 | Pass | User |
| ... | ... | ... | ... | ... | ... | 06/19/2012 11:25 am | ... | 1 | Pass | User |
| ... | ... | ... | ... | ... | ... | 06/19/2012 11:25 am | ... | 1 | Pass | User |
| ... | ... | ... | ... | ... | ... | 06/19/2012 11:24 am | ... | 1 | Pass | User |
| ... | ... | ... | ... | ... | ... | 06/19/2012 11:24 am | ... | 1 | Pass | User |
| ... | ... | ... | ... | ... | ... | 06/19/2012 11:24 am | ... | 1 | Pass | User |
| ... | ... | ... | ... | ... | ... | 06/19/2012 11:24 am | ... | 1 | Pass | User |

SMART Link®

The wireless gateway, located at each site, allows the SMART Guard® and SMART Shield® to communicate to the SMART Command Center® Website via Ethernet, WiFi, or Cellular Communication.



SMART Shield®

A wireless handheld thermocouple used to take food temperatures throughout the flow of food: receiving, preparation, cooking, holding, serving, and cooling.



- Automatically downloads food menu items daily from the SMART Command Center® Website.
- Confirms and documents critical temperature limits and prompts immediate corrective actions when these limits fall below, or exceed, the acceptable temperature range.

SMART Shield® Pro

A touchscreen Android™ device with a custom SMART Shield® Pro Application, and an option for a Bluetooth wireless temperature probe.



- All the capabilities of our SMART Shield® with automatic downloading of menu items, user information, task lists, and inspection forms via our SMART Command Center® Website.
- This unit may be taken from site to site and can be used as a verification device and inspection tool.



All products provide, real time temperature documentation that **exceeds all SDOH, TJC, CDC, CMS, FDA, and HACCP Reporting Compliance** standards via a secure wireless network to the **SMART Command Center®** Website. The entire SMART Temps Temperature Management System exceeds Joint Commission Standards and helps your team put patient safety at the forefront with a real return on investment.

Contact SMART Temps today for your 30-day trial.

877-272-3111 | SMART-Temps.com